

NAME OF PROTOCOL

The tissue clearing protocol for organoid (Cat. No.SHOC-001)

CONDITION OF SAMPLE

Length 1 mm x Width 1 mm x Height 1 mm

CODE OF PROTOCOL

C1003

REVISION OF PROTOCOL

1.1.2 (2019.08.21)

Preparation | Planning your test.



When we designed the protocol, we considered not only the effect of the clearing but also the working time of the researchers. *Enjoy the tissue clearing!*

Preparation | Takes the solutions.

A-1. All of the solutions might be stored at 4°C.

A-2. Crystallization or precipitation while solutions are stored at 4°C; It is no effect for the clearing process. It will be liquid again after incubating the solutions at 37°C (Refer to preparation step B-3).

A-3. Do not use the individual solution from the other kit. Even if the names of solutions are the same, the component compositions are not the same. Each solution has a unique component composition depend on the purpose of the kit.

1	Starting Solution	30 ml
2	Organoid Clearing Solution	30 ml
+	Mounting & Storage Solution	60 ml

- The Mounting & Storage Solution (Cat. No. SHMS-060) is **not included in the kit**. Only the Starter's kits contain the Mounting & Storage Solution.
- The solutions may be **crystallized or precipitated**, incubate it at 37°C for 1~2 h before use.

Preparation | Fixing the sample (organoid tissue).

B-1. Incubate the sample with 4% PFA at 4°C for 15 min.

B-2. Wash the sample with 1 X PBS in shaking incubator at room temperature for 10 min X 3 times.

B-3. Incubate the **2 Organoid Clearing Solution** and **+ Mounting & Storage Solution** at 37°C for 1~2 h before use.

Protocol | Clearing the fixed sample.

1. Incubate the sample with **1 500 µl Starting Solution** at 4°C for **until the sample sinks (at least overnight or more)**.

2. Incubate the sample with **2 500 µl Organoid Clearing Solution** in shaking incubator at 50 rpm/37°C for 6 h.

3. Wash the sample with **distilled water** in shaking incubator at 50 rpm/4°C for 10 min X 3 times.

⚠ The sample may be **back to be opaque reversely and swell**.

It is no effect for the clearing process, the sample will be cleared again at the end of protocol.

4. Incubate the sample with **+ 500 µl Mounting & Storage Solution** in shaking incubator at 50 rpm/37°C

for at least 4 h or more.

Clearing Tips | Troubleshootings

- C-1. If the sample contains air bubbles ➡ Centrifuge the sample at 3,500 rpm/24°C for 30 min.
- C-2. If the sample is not entirely clearing ➡ Repeat step #3 and then from #2 to step #4.
- C-3. If the rpm is not specified ➡ Usually operate the shaking incubator gently.
- C-4. **Never wash the sample with PBS** instead of distilled water at step #3.
- C-5. Use the vial or tube for tissue clearing better than the slide chamber.

Storage & Imaging Tips

- D-1. Store the cleared sample in **+ Mounting & Storage Solution** at the room temperature (20~25°C).
- D-2. **Take images within 7 days after the clearing** for the best result.
- D-3. Take images on the microscope. We recommend using the Light Sheet Fluorescence Microscope (LSFM) or Confocal Laser Scanning Microscope (CLSM). And analyze and visualize the images with the microscopy image analysis software.
- D-4. **+ Mounting & Storage Solution** the **solvent-free** material that is safe for using the Light Sheet Fluorescence Microscope (LSFM).
- D-5. **Refractive Index (RI)** of **+ Mounting & Storage Solution** is **1.45**.
- D-6. **Be careful of making bubbles** while filling the microscope chamber with the sample and **+ Mounting & Storage Solution**. The bubbles may disturb taking imaging.

End of Protocol

Contact Us | Technical support

Binaree, Inc. (Headquarters)
47 Gyeongdaero17-gil Buk-gu,
STE#608 IT Convergence Bldg(115),
Daegu, 41566, Republic of Korea.

◦ Website : binaree.com ◦ Tel. +82-(0)53-291-5012
◦ Email : lab@binaree.com ◦ Fax. +82-(0)53-382-5012